




# Study the dynamics of drug delivery

By recording a sequence of images and letting the system automatically do the analysis for you


Step 1. Drug application by iontophoresis.



Step 2. Capture a sequence of images through the ring-shape iontophoresis electrode and watch the result.



Step 3. Display the time-sequence of average TiVi-values.



Step 4. Generate a dose-response curve.

*WheelsBridge*

PIONEERS IN TISSUE VIABILITY IMAGING



WheelsBridge AB, Lövsbergsvägen 13, S-589 37 Linköping, Sweden  
Phone: +46-708-765-190, E-mail: info@wheelsbridge.se, WEB-site: www.wheelsbridge.se