


# Skin Colour Tracker


Displays alterations in local skin appearance by tracking its colour changes




Step 1. Display a photo and select a Region of Interest (ROI).



Step 2. Select Colours of Interest (COI).




Step 3. Watch how these colours are represented in a sequence of photos.



Step 4. Display the trend.

*WheelsBridge*

PIONEERS IN TISSUE VIABILITY IMAGING



WheelsBridge AB, Lövsbergsvägen 13, S-589 37 Linköping, Sweden  
 Phone: +46-708-765-190, E-mail: info@wheelsbridge.se, WEB-site: www.wheelsbridge.se