


# Let the Express Wizard Guide You

through the process of analyzing your data captured by the TiVi600 tissue Viability Imager

Step 1. Apply the test substance and capture a sequence of images.



Step 2. Load the images for analysis.



Step 3. Set threshold parameters in the Express Wizard.



Step 4. View the results in the Erythema Intensity/Area windows.



Using the Express Wizard significantly speeds up the analysis process and makes your work more productive.

*WheelsBridge*

PIONEERS IN TISSUE VIABILITY IMAGING



WheelsBridge AB, Lövsbergsvägen 13, S-589 37 Linköping, Sweden  
Phone: +46-708-765-190, E-mail: info@wheelsbridge.se, WEB-site: www.wheelsbridge.se